# Destination®

Carotid Guiding Sheath / Introduceur carotidien / Halsschlagader-  
Führungsschleuse / Vaina guía carotídea / Bainha-guia Carótida /  
Introduttore carotideo / Carotisgeleidesheath / Styrinroducer för  
halspulsådern / Carotisledehylster / Καρωτιδική κατευθυντήρια θήκη

Read These Instructions Before Use / Lire ces instructions avant utilisation /  
Bitte diese Anleitungen vor dem Gebrauch lesen / Leer estas instrucciones  
antes del uso / Leia estas Instruções antes de Utilizar / Leggere le presenti  
istruzioni prima dell'uso / Lees deze aanwijzingen vóór gebruik /  
Läs dessa anvisningar före användning / Læs denne brugsanvisning inden  
anvendelse / Διαβάστε αυτές τις οδηγίες πριν από τη χρήση

Instructions For Use .....	3
Mode d'emploi.....	8
Gebrauchsanleitung .....	14
Instrucciones de uso.....	20
Instruções para Utilização .....	26
Istruzioni per l'uso .....	32
Gebruiksaanwijzing .....	38
Bruksanvisning .....	44
Brugsanvisning .....	50
Οδηγίες χρήσης.....	56

 **TERUMO** CE0197

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**TERUMO MEDICAL CORPORATION**  
950 ELKTON BLVD, ELKTON, MD 21921 USA  
MADE IN USA

**EC REP**

**TERUMO EUROPE N.V.**  
INTERLEUVENLAAN 40, 3001 LEUVEN, BELGIUM

**TERUMO CORPORATION**  
TOKYO 151-0072, JAPAN

**AU Sponsor**

**TERUMO AUSTRALIA PTY LTD**  
MACQUARIE PARK NSW 2113 AUSTRALIA

## DESCRIPTION

The Destination<sup>®</sup> Carotid Guiding Sheath is designed to perform as a guiding catheter and an introducer sheath. The sheath is coil reinforced, has a radiopaque marker and approximately 15 cm (150 mm) of the distal tip is hydrophilically coated. It comes packaged with a dilator and choice of hemostatic valve types.

## INTENDED USE


The Destination Carotid Guiding Sheath is designed to be used for the introduction of interventional and diagnostic devices into the human vasculature, including but not limited to the carotid arteries.

## TARGET PATIENT POPULATION

Destination Guiding Sheath is for adult population undergoing diagnostic and/or interventional endovascular procedures. This device has not been tested on vulnerable patient populations (e.g., Pregnant women).

## COMPONENT DESCRIPTION


Refer to product labeling for appropriate system components. All components that enter blood vessels are radiopaque.

The contents of this package may differ from that shown in the figure above.

## Sheath

This highly flexible, reinforced sheath is designed to resist kinking. It incorporates a radiopaque marker located approximately 5 mm proximal to the tip (Fig. 1). The Destination Carotid Guiding Sheath is available in multiple distal segments shapes. The sheath has an atraumatic tip which is hydrophilically coated. Coating length is approximately 15 cm (150 mm) as measured from the distal end.


## Dilator

The precise fit of the dilator in the sheath allows for simultaneous motion of both dilator and sheath. The radiopacity of the dilator has been greatly enhanced to increase the visibility. The distal segment of the dilator has been optimized to extend approximately 5 cm (50 mm) beyond the tip of the sheath.


## Valve

Each Destination Carotid Guiding Sheath comes packaged with either a Cross Cut Valve assembly (CCV) or a Y-shaped Tuohy-Borst Valve (TBV). Both valve types are equipped with a three way stopcock connected by a side tube.

## Dilator Retaining Clip

The clip provides a positive mechanism for locking the dilator hub and the TBV together. The clip prevents axial movement of the dilator while the sheath is being advanced (Fig. 2).

## CONTRAINDICATIONS

There are no contraindications for the Destination Guiding Sheath.

## PRECAUTIONS

- This device must only be used by a trained physician.
- This device is for single use only. Do not re-sterilize or re-use.
- Destination is a single-use device. Harms that may result from reuse include, but are not limited to:
  - I. Iatrogenic infection
  - II. Foreign body reaction
  - III. Allergic/Toxic reaction
- Contents are sterile, non-toxic and non-pyrogenic in an unopened undamaged package.
- Do not use if the package or product is stained or damaged. Use the device immediately after opening the package and dispose of the device after use.
- Before use make sure the sheath size (Fr.) is appropriate for the access vessel and the interventional/ diagnostic device to be used.
- The entire procedure, from skin incision to sheath removal, must be carried out aseptically.
- Do not use a power injector through the side tube and 3-way stopcock. Excessive leakage may occur through the CCV valve with high/rapid flow injections such as an injection of contrast to provide a comprehensive image of the aortic arch.
- When inflating a balloon at, or close to, the sheath tip, be sure not to inflate it inside the distal end of the sheath (Fig. 6). The radiopaque marker is located approx. 5 mm proximal to the tip and therefore does not actually mark the true distal tip of the sheath (Fig. 1).
- Do not heat or bend the sheath tip. Damage to the sheath may result.

### **CAUTION**

- Federal (U.S.A) law restricts this device to sale by or on the order of a physician.


## INSTRUCTIONS FOR USE

Instructions for using the Tuohy-Borst Valve (TBV) are covered in sections I and III.

Instructions for using the Cross Cut Valve (CCV) are covered in sections II and III.

### I. Using the Tuohy-Borst Valve (TBV)

1. Ensure that the valve and the sheath are tightly connected.
2. Connect a flushing line to the 3-way stopcock of the guiding sheath. Fill the sheath assembly completely with heparinized saline, removing all air.
3. Wet the distal 15 cm (150 mm) of the sheath with heparinized saline. To maintain lubricity, this surface must be kept completely wet.
4. Prime the dilator using a syringe with heparinized saline.
5. Insert the dilator completely into the sheath.
6. Tighten the Tuohy-Borst Valve around the dilator by turning the cap clockwise. If desired, clip the dilator hub and the TBV together using the dilator retaining clip provided (Figs. 2 & 3). Be sure to note the correct orientation of the dilator retaining clip to the dilator hub and TBV.
7. To “unlock” the dilator from the TBV and sheath, remove the clip and loosen the TBV by turning the cap counter clockwise.

#### CAUTION

- The Tuohy-Borst Valve (TBV) must be tightly closed to minimize the risk of hemorrhage or air embolism.

### II. Using the Cross Cut Valve (CCV)



1. Ensure that the valve and the sheath are tightly connected.
2. Connect a flushing line to the 3-way stopcock of the guiding sheath. Fill the sheath assembly completely with heparinized saline, removing all air.
3. Wet the distal 15 cm (150 mm) of the sheath with heparinized saline. To maintain lubricity, this surface must be kept completely wet.
4. Prime the dilator using a syringe with heparinized saline.
5. Insert the dilator completely into the sheath.

#### CAUTION

- Insert the dilator into the center of the sheath valve. Forced insertion of the dilator which misses the center of the sheath valve may cause damage and result in blood leakage.
6. The CCV connects with the male hub of the dilator, and locks in place by means of grip. Snap the dilator hub into the top of the CCV hub to lock them together (Fig. 4).
  7. To “unlock” the dilator from the CCV and sheath, bend the dilator hub downward (Fig. 5).

#### CAUTION

- Be sure that the dilator is held firmly in place in the sheath and CCV. If the dilator is not secure, only the sheath will advance into the vessel and the tip of the sheath may damage the vessel. Advancing the sheath alone may damage the vessel.

### III. For All Valve Types (TBV and CCV)

1. Make a small skin incision at the puncture site with a surgical knife.
2. Insert a cannula into the vessel.
3. Insert a guidewire through the cannula into the vessel.

#### **CAUTION**

- Advance or withdraw the guidewire slowly. If resistance is met, do not advance or withdraw the guidewire until the cause of resistance has been determined.
4. Remove the cannula over the guidewire.
  5. Insert the dilator and sheath together over the guidewire and into the blood vessel and advance it to the target vessel.

#### **CAUTION**

- Advance or withdraw the guidewire slowly. If resistance is met, do not advance or withdraw the guidewire until the cause of resistance has been determined.
6. Slowly remove the dilator leaving the sheath in the vessel. If injection or sampling is necessary at this point, remove the guidewire only and use the dilator hub as an injection port before removing it.

#### **CAUTION**


- TBV Type: Loosen the TBV before removing the dilator. Then tightly close the TBV to minimize the risk of hemorrhage or air embolism.
  - CCV Type: Be sure to remove the dilator from the sheath slowly. Rapid withdrawal of the dilator may result in the incomplete closing of the CCV valve, resulting in blood flow through the valve. If this occurs, replace the dilator into the sheath and remove again slowly.
7. The radiopaque marker will identify the sheath tip location under fluoroscopy. The radiopaque marker is located approximately 5 mm proximal to the sheath's tip (fig.1).
  8. Insert an interventional/diagnostic device through the sheath and into the blood vessel, then advance it to the desired location.

#### **CAUTION**

- Advance or withdraw the sheath slowly. If resistance is met, do not advance or withdraw the sheath until the cause of resistance has been determined.
- Before removing or inserting the interventional/ diagnostic device through the sheath, aspirate blood from the 3-way stopcock to remove any fibrin deposition which may have accumulated in or on the tip of the sheath.
- When puncturing, suturing or incising the tissue near the sheath, be careful not to damage the sheath. Do not put a clamp on the sheath or bind it with a thread.
- When inflating a balloon at, or close to, the sheath tip, be sure not to inflate it inside the distal end of the sheath (Fig. 6). The

radiopaque marker is located approximately 5 mm proximal to the tip (Fig. 1) and therefore does not actually mark the true distal tip of the sheath.

9. The side tube may also be used as a continuous infusion site by connecting an infusion line to the 3-way stopcock.
10. When inserting, manipulating or withdrawing an interventional/ diagnostic device from the sheath, always be sure to hold the sheath in place.

### **CAUTION**

- Do not suture the sheath tubing since this may restrict access / flow through the sheath.
  - Do not use a power injector through the side tube and 3-way stopcock. Excessive leakage may occur through the CCV valve with high/rapid flow injections such as an injection of contrast to provide a comprehensive image of the aortic arch.
11. After the intended procedure is completed, remove all of the devices.

### **STORAGE**

Do not store at extreme temperatures and humidity. Avoid direct sunlight.

### **POTENTIAL ADVERSE EVENTS AND COMPLICATIONS**

Adverse events that may occur:

- Allergic/Toxic reactions.
- Arterial dissection.
- Arterial Perforation.
- Blood Loss (Hemorrhage/Hematoma).
- Damage to tissue at entry site.
- Dosage error, PRN (as needed) drug leaks out as it is injected in patient.
- Embolism.
- Embolism (Foreign Matter-product material).
- Embolism, Air.
- Foreign body reaction (also formation and/or introduction of thrombus into vasculature).
- Foreign body reaction.
- Harm from incorrect deployment of interventional device.
- Iatrogenic infection – Systemic (circulating blood).
- Iatrogenic infection – Local (access site).
- Kidney/Renal insufficiency, Acute (excess use of contrast/media due to prolonged procedure as a result of device damage).
- Procedural delay – replacement of existing device.
- Procedure Delay (out of the box failure) - Defect easily identifiable. Device not used. Delay in setup.
- Public or Environmental health hazard – cross contamination to health care worker.
- Vascular Spasm.
- Vascular Surgical Intervention.
- Cardiac Tamponade
- Compartment syndrome
- Cut, Abrasion or Puncture wound
- Evulsion
- Excessive Exposure to fluoroscopy
- Kidney/ Renal Failure

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**TERUMO** and **TERUMO** are registered in the U.S. Patent and Trademark Office by **TERUMO CORPORATION**. Destination is a trademark of **TERUMO MEDICAL CORPORATION**.
**TERUMO** und **TERUMO** wurden im US Patent and Trademark Office von **TERUMO CORPORATION** registriert. Destination wurde im US Patent & Trademark Office von **TERUMO MEDICAL CORPORATION** registriert.
**TERUMO** e **TERUMO** estão registradas no U.S. Patent and Trademark Office (Escritório para Registo de Marcas e Patentes dos Estados Unidos) pela **TERUMO CORPORATION**. Destination é uma marca comercial da **TERUMO MEDICAL CORPORATION**.
**TERUMO** en **TERUMO** zijn door **TERUMO CORPORATION** gedeponeerd bij het U.S. Patent and Trademark Office. Destination is een handelsmerk van **TERUMO MEDICAL CORPORATION**.