

# ***FineCross® M3***

Coronary Micro-Guide Catheter



## OPTIMAL GUIDEWIRE SUPPORT AND SUPERIOR **CROSSABILITY**<sup>1</sup>

**FINECROSS® M3 Microcatheter is designed to provide the ultimate balance of properties for access and crossing of complex lesions.<sup>1</sup>**

- **1.7 Fr tapered distal tip** is designed for improved crossability and guidewire handling<sup>1</sup>
- **Next generation hydrophilic coating**, with its unmatched and most durable lubricity, provides excellent trackability and crossability<sup>1</sup>
- **15 cm floppy distal segment** allows for improved trackability
- **Tapered stainless steel braid construction** is designed for optimal guidewire support

# PUSHING BOUNDARIES

Terumo Interventional Systems is **committed to your success** with innovative procedural solutions and ongoing support for your most challenging cases.


We are relentlessly seeking new ways to help you achieve **better outcomes for more patients.**

## ORDERING INFORMATION

### FINECROSS M3 Coronary Micro-Guide Catheter

PRODUCT CODE	EFFECTIVE CATHETER LENGTH (cm)	DISTAL ENTRY TIP PROFILE (Fr/mm)	DISTAL OUTER DIAMETER (Fr/mm)	PROXIMAL OUTER DIAMETER (Fr/mm)	DISTAL HYDROPHILIC COATING (cm)	DISTAL INNER DIAMETER (mm)	PROXIMAL INNER DIAMETER (in/mm)	MAXIMUM GUIDEWIRE SIZE (in/mm)
<b>M3-1430</b>	130	1.7/0.55	1.8/0.61	2.5/0.84	75	0.45	0.021/0.55	0.014/0.36
<b>M3-1450</b>	150	1.7/0.55	1.8/0.61	2.5/0.84	95	0.45	0.021/0.55	0.014/0.36


**Note:** Product specifications for FINECROSS MG, codes 35-1430 and 35-1450, are 1.8/0.60 (Fr/mm) for Distal Entry Tip Profile, 2.6/0.87 (Fr/mm) for Proximal Outer Diameter, and 70 and 90 (cm) for Distal Hydrophilic Coating respectively.



## SUPPORT

### for Side-branch or Acute Bifurcation Lesions



Use FINECROSS M3 Microcatheter with RUNTHROUGH® NS Extra Floppy Coronary Guidewire



## CROSSABILITY

### for Distal, Tortuous and Tight Stenotic Lesions

Use FINECROSS M3 Microcatheter with RUNTHROUGH® NS HYPERCOAT™ Coronary Guidewire

**FIND OUT MORE**  Phone: 800.888.3786  [terumo.com](https://www.terumo.com)

**RX ONLY.** Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

**References:**  
1. Data on File.

©2023 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of their respective owners. PM-06948

**TERUMO**  
INTERVENTIONAL  
SYSTEMS